

ATF3 Rabbit mAb


ABZ3261

Product Overview

Name	ATF3 Rabbit mAb
Accession(Primary)	P18847
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:200; ICC/IF-1:50-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (P18847)

Gene Name	ATF3
About Gene Symbol	



Western blot analysis of ATF3 in MCF-7, U251, HL-60 lysates using ATF3 antibody.