

ATF6 Rabbit mAb


ABZ3262

Product Overview


Name	ATF6 Rabbit mAb
Accession(Primary)	P18850
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 75 kDa; Observed MW: 90-100 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20

Target information (P18850)

Gene Name	ATF6
About Gene Symbol	



Western blot analysis of ATF6 in A549, HL-60, U2OS lysates using ATF6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using ATF6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.