

BubR1 Rabbit mAb


ABZ3264

Product Overview


Name	BubR1 Rabbit mAb
Accession(Primary)	O60566
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 120 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:200-1:500; IP-1:20-1:50

Target information (O60566)

Gene Name	BUB1B
About Gene Symbol	BUB1B (BUB1 Mitotic Checkpoint Serine/Threonine Kinase B) is a serine/threonine protein kinase involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of BubR1 in HeLa, A549, MCF-7, U251 lysates using BUBR1 antibody.



Immunocytochemistry analysis of BubR1 (green) in HeLa using BubR1 antibody, and DAPI (blue)