

BubR1 Rabbit mAb


ABZ3265

Product Overview


Name	BubR1 Rabbit mAb
Accession(Primary)	O60566
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 120 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (O60566)

Gene Name	BUB1B
About Gene Symbol	BUB1B (BUB1 Mitotic Checkpoint Serine/Threonine Kinase B) is a serine/threonine protein kinase involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of BubR1 in CHO-K1, Ramos lysates using BUBR1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using BubR1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.