


Cadherin 10 Rabbit mAb


ABZ3266

Product Overview

Name	Cadherin 10 Rabbit mAb
Accession(Primary)	Q9Y6N8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (Q9Y6N8)

Gene Name CDH10**About Gene Symbol**



Western blot analysis of cadherin 10 in U2OS, C6
lysates using cadherin 10 antibody.