

Calumenin Rabbit mAb


ABZ3268

Product Overview


Name	Calumenin Rabbit mAb
Accession(Primary)	O43852
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:10-1:100

Target information (O43852)

Gene Name	CALU
About Gene Symbol	



Western blot analysis of Calumenin in HeLa, A549, C6 lysates using Calumenin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Calumenin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.