

CBX4 Rabbit mAb


ABZ3270

Product Overview


Name	CBX4 Rabbit mAb
Accession(Primary)	O00257
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:200-1:500; IP-1:20-1:50

Target information (O00257)

Gene Name	CBX4
About Gene Symbol	



Western blot analysis of CBX4 in 3T3 lysates using CBX4 antibody.



Western blot analysis of CBX4 in HeLa, A549, HL-60, U2OS, C6 lysates using CBX4 antibody.