

CD11a Rabbit mAb

ABZ3271


Product Overview

Name	CD11a Rabbit mAb
Accession(Primary)	P24063
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 128 kDa; Observed MW: 160 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:500-1:2000


Target information (P24063)

Gene Name Itgal

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded mouse lung using Integrin alpha L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD11a in A549, HL-60
lysates using CD11a antibody.