

CD229 Rabbit mAb

ABZ3272


Product Overview

Name	CD229 Rabbit mAb
Accession(Primary)	Q9HBG7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500

Target information (Q9HBG7)

Gene Name LY9

About Gene Symbol



Western blot analysis of CD229 in HeLa, A549, HL-60, U251 lysates using CD229 antibody.

