

CEP55 Rabbit mAb


ABZ3275

Product Overview


Name	CEP55 Rabbit mAb
Accession(Primary)	Q53EZ4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:100-1:200

Target information (Q53EZ4)

Gene Name	CEP55
About Gene Symbol	CEP55 (Centrosomal Protein 55) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of CEP55 (green) in HeLa using CEP55 antibody, and DAPI (blue).



Western blot analysis of CEP55 in HeLa, A549, HL-60, U2OS, C6 lysates using CEP55 antibody.