

CLIC1 Rabbit mAb


ABZ3277

Product Overview


Name	CLIC1 Rabbit mAb
Accession(Primary)	O00299
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (O00299)


Gene Name	CLIC1
About Gene Symbol	



Immunocytochemistry analysis of CLIC1 (green) in Hela using CLIC1 antibody, and DAPI (blue).



Western blot analysis of CLIC1 in Raw264.7, mouse kidney, mouse spleen lysates using CLIC1 antibody.



Western blot analysis of CLIC1 in HeLa, CHO-K1, C6, Ramos lysates using CLIC1 antibody