

Cytochrome P450 2D6 Rabbit mAb

ABZ3281


Product Overview

Name	Cytochrome P450 2D6 Rabbit mAb
Accession(Primary)	P10635
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 56 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P10635)

Gene Name CYP2D6


About Gene Symbol



Western blot analysis of CYP2D6 in HeLa, A549, HL-60, U251, U87-MG lysates using Cytochrome P450 2D6 antibody.



Immunohistochemistry analysis of paraffin-embedded liver cancer using CYP2D6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of CYP2D6 (green) in HeLa using CYP2D6 antibody, and DAPI (blue)