


# Cystatin B Rabbit mAb


ABZ3282

## Product Overview


|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Cystatin B Rabbit mAb                      |
| <b>Accession(Primary)</b>   | P04080                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Rat                                 |
| <b>Calculated MW</b>        | Calculated MW: 11 kDa; Observed MW: 11 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200        |

## Target information (P04080)

**Gene Name** CSTB**About Gene Symbol**



Western blot analysis of CystatinB in HeLa, A549, HL-60 lysates using CystatinB antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CystatinB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.