

DGKZ/DGK-zeta Rabbit mAb


ABZ3283

Product Overview

Name	DGKZ/DGK-zeta Rabbit mAb
Accession(Primary)	Q13574
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 124,114 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:20-1:50; FC-1:200-1:500

Target information (Q13574)

Gene Name	DGKZ
About Gene Symbol	



Western blot analysis of DGKZ/DGKzeta in HeLa, A549, HL-60, U2OS lysates using DGKZ antibody.