

DGKZ/DGK-zeta Rabbit mAb


ABZ3284

Product Overview


Name	DGKZ/DGK-zeta Rabbit mAb
Accession(Primary)	Q13574
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 124,114 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; FC-1:200-1:500; IP-1:20-1:50

Target information (Q13574)

Gene Name	DGKZ
About Gene Symbol	



Western blot analysis of DGKZ/DGKzeta in mouse brain, rat brain lysates using DGKZ antibody.



Western blot analysis of DGKZ/DGKzeta in HeLa, A549, HL-60, U2OS, U87-MG lysates using DGKZ/DGKzeta antibody.