


DLX3 Rabbit mAb

ABZ3285

Product Overview

Name	DLX3 Rabbit mAb
Accession(Primary)	O60479
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O60479)

Gene Name DLX3**About Gene Symbol**

Western blot analysis of DLX3 in A549, U2OS
lysates using DLX3 antibody.

