

EBP1 Rabbit mAb


ABZ3287

Product Overview


Name	EBP1 Rabbit mAb
Accession(Primary)	Q9UQ80
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (Q9UQ80)


Gene Name	PA2G4
About Gene Symbol	PA2G4 (Proliferation-Associated 2G4) is a membrane receptor involved in cell cycle regulation. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify th



Western blot analysis of EBP1 in 3T3, mouse liver lysates using EBP1 antibody.



Western blot analysis of EBP1 in CHO-K1, C6, Ramos lysates using EBP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using EBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.