

eIF4B Rabbit mAb


ABZ3289

Product Overview


Name	eIF4B Rabbit mAb
Accession(Primary)	P23588
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 69 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:100-1:200

Target information (P23588)

Gene Name	EIF4B
About Gene Symbol	eIF4B (Eukaryotic Translation Initiation Factor 4B) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of eIF4B in HeLa, A549, HL-60, U2OS lysates using eIF4B antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using eIF4B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.