

Entactin Rabbit mAb

ABZ3290


Product Overview

Name	Entactin Rabbit mAb
Accession(Primary)	P14543
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 136 kDa; Observed MW: 136 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:1000


Target information (P14543)

Gene Name NID1

About Gene Symbol



Western blot analysis of Entactin/NID in HeLa, A549, HL-60, U251, U87-MG lysates using Entactin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Entactin/NID antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.