

ErbB3 / HER3 Rabbit mAb


ABZ3291

Product Overview

Name	ErbB3 / HER3 Rabbit mAb
Accession(Primary)	P21860
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 148 kDa; Observed MW: 185 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:20-1:50

Target information (P21860)

Gene Name	ERBB3
About Gene Symbol	



Western blot analysis of ErbB 3 in A549, MCF-7
lysates using ErbB 3 antibody.