

EYA2 Rabbit mAb

ABZ3293


Product Overview

Name	EYA2 Rabbit mAb
Accession(Primary)	O00167
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O00167)

Gene Name EYA2

About Gene Symbol



Western blot analysis of EYA2 in HeLa, A549, U2OS lysates using EYA2 antibody.

