

FADS1 Rabbit mAb


ABZ3295

Product Overview



Name	FADS1 Rabbit mAb
Accession(Primary)	O60427
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:100-1:1000

Target information (O60427)


Gene Name	FADS1
About Gene Symbol	




Immunocytochemistry analysis of FADS1 (green) in K562 using FADS1 antibody, and DAPI (blue).



Western blot analysis of FADS1 in rat brain, rat liver lysates using FADS1 antibody.



Western blot analysis of FADS1 in rat Brain
lysates using FADS1 antibody



Immunohistochemistry analysis of paraffin-
embedded Human lung cancer using FADS1
antibody. High-pressure and temperature Sodium
Citrate pH 6.0 was used for antigen retrieval.