

GALT Rabbit mAb

ABZ3301


Product Overview

Name	GALT Rabbit mAb
Accession(Primary)	P07902
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100


Target information (P07902)

Gene Name

GALT

About Gene Symbol

Western blot analysis of GALT in mouse liver, mouse kidney lysates using GALT antibody.



Western blot analysis of GALT in HeLa, C6, rat Brain, Ramos lysates using GALT antibody.