

GMEB1 Rabbit mAb


ABZ3305

Product Overview

Name	GMEB1 Rabbit mAb
Accession(Primary)	Q9Y692
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:20-1:50; FC-1:200-1:500; IP-1:20-1:50

Target information (Q9Y692)

Gene Name	GMEB1
About Gene Symbol	



Western blot analysis of GMEB1 in A549, HL-60, U2OS lysates using GMEB1 antibody.