

GNAI2 Rabbit mAb


ABZ3306

Product Overview


Name	GNAI2 Rabbit mAb
Accession(Primary)	P04899
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; FC-1:10-1:100; ICC/IF-1:50-1:100

Target information (P04899)


Gene Name	GNAI2
About Gene Symbol	




Immunocytochemistry analysis of GNAI2 (green)
in HeLa using GNAI2 antibody, and DAPI (blue).



Western blot analysis of GNAI2 in U87-MG lysates
using GNAI2 antibody.



Western blot analysis of GNAI2 in C6, mouse Brain lysates using GNAI2 antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using GNAI2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.