


HCK Rabbit mAb

ABZ3308


Product Overview

Name	HCK Rabbit mAb
Accession(Primary)	P08631
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:2000-1:5000; ICC/IF-1:50-1:100

Target information (P08631)

Gene Name HCK**About Gene Symbol**

Western blot analysis of Hck in HeLa, A549, HL-60, U2OS lysates using Tyrosine Protein Kinase HCK antibody.



Immunocytochemistry analysis of Hck (green) in
hela using Hck antibody, and DAPI(blue)