

hHR23A Rabbit mAb


ABZ3309

Product Overview


Name	hHR23A Rabbit mAb
Accession(Primary)	P54725
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 40 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P54725)

Gene Name	RAD23A
About Gene Symbol	RAD23A (RAD23 Homolog A, Nucleotide Excision Repair Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of hHR23A in 3T3 lysates using hHR23A antibody.



Western blot analysis of hHR23A in A549, HL-60, U2OS lysates using hHR23A antibody.