

## HINT1 Rabbit mAb


ABZ3310

### Product Overview


|                             |                                                             |
|-----------------------------|-------------------------------------------------------------|
| <b>Name</b>                 | HINT1 Rabbit mAb                                            |
| <b>Accession(Primary)</b>   | P49773                                                      |
| <b>Host Species</b>         | Rabbit                                                      |
| <b>Clonality</b>            | Monoclonal                                                  |
| <b>Isotype</b>              | IgG                                                         |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                           |
| <b>Calculated MW</b>        | Calculated MW: 14 kDa; Observed MW: 14 kDa                  |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:500-1:1000 |

### Target information (P49773)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | HINT1 |
| <b>About Gene Symbol</b> |       |



Immunocytochemistry analysis of HINT1 (green) in HeLa using HINT1 antibody, and DAPI (blue).



Western blot analysis of HINT1 in A549, HL-60, U251 lysates using HINT1 antibody.