

Histone H2A.X Rabbit mAb


ABZ3311

Product Overview


Name	Histone H2A.X Rabbit mAb
Accession(Primary)	P16104
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20-1:50

Target information (P16104)


Gene Name	H2AX
About Gene Symbol	



Immunocytochemistry analysis of Histone H2A.X (green) in HeLa using Histone H2A.X antibody, and DAPI (blue).



Western blot analysis of Histone H2A.X in HeLa, Ramos lysates using Histone H2A.X antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using Histone H2A.X antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.