

HMG4 Rabbit mAb


ABZ3312

Product Overview


Name	HMG4 Rabbit mAb
Accession(Primary)	O15347
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20

Target information (O15347)


Gene Name	HMGB3
About Gene Symbol	



Western blot analysis of HMG4 in 3T3 lysates using HMG4 antibody.



Western blot analysis of HMG4 in HeLa, A549, HL-60, U2OS, C6 lysates using HMG4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HMG4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.