

ICAM3 Rabbit mAb


ABZ3314

Product Overview


Name	ICAM3 Rabbit mAb
Accession(Primary)	P32942
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 110-160 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200

Target information (P32942)

Gene Name	ICAM3
About Gene Symbol	



Western blot analysis of ICAM3 in HeLa, A549, HL-60, U2OS, C6 lysates using ICAM3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ICAM3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.