

IGF1 Receptor Rabbit mAb


ABZ3315

Product Overview

Name	IGF1 Receptor Rabbit mAb
Accession(Primary)	P08069
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 155 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (P08069)

Gene Name	IGF1R
About Gene Symbol	



Western blot analysis of IGF1 Receptor in HeLa, A549, MCF-7, U251 lysates using IGF1 Receptor antibody.