


# IGF1 Receptor Rabbit mAb


ABZ3316

## Product Overview

|                             |                                             |
|-----------------------------|---------------------------------------------|
| <b>Name</b>                 | IGF1 Receptor Rabbit mAb                    |
| <b>Accession(Primary)</b>   | P08069                                      |
| <b>Host Species</b>         | Rabbit                                      |
| <b>Clonality</b>            | Monoclonal                                  |
| <b>Isotype</b>              | IgG                                         |
| <b>Species Reactivity</b>   | Human                                       |
| <b>Calculated MW</b>        | Calculated MW: 155 kDa; Observed MW: 95 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:200-1:500            |

## Target information (P08069)

**Gene Name** IGF1R**About Gene Symbol**



Immunohistochemistry analysis of Human lung cancer using IGF1 Receptor antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.