


IKZF3 Rabbit mAb

ABZ3318


Product Overview

Name	IKZF3 Rabbit mAb
Accession(Primary)	Q9UKT9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (Q9UKT9)

Gene Name IKZF3**About Gene Symbol**

Western blot analysis of IKZF3 in HL-60, U251 lysates using IKZF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IKZF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.