

IL-26 Rabbit mAb


ABZ3319

Product Overview


Name	IL-26 Rabbit mAb
Accession(Primary)	Q9NPH9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20; IF-1:50-1:200

Target information (Q9NPH9)

Gene Name	IL26
About Gene Symbol	



Western blot analysis of IL-26 in HeLa, A549, HL-60, U251 lysates using IL-26 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IL-26 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.