


KIFAP3 Rabbit mAb

ABZ3323


Product Overview

Name	KIFAP3 Rabbit mAb
Accession(Primary)	Q92845
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 91 kDa; Observed MW: 91 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (Q92845)

Gene Name KIFAP3**About Gene Symbol**

Western blot analysis of KIFAP3 in A549, HL-60, C6 lysates using KIFAP3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using KIFAP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.