

Ku70 Rabbit mAb


ABZ3325

Product Overview


Name	Ku70 Rabbit mAb
Accession(Primary)	P12956
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (P12956)


Gene Name	XRCC6
About Gene Symbol	XRCC6 (X-Ray Repair Cross Complementing 6) is an ATP-dependent DNA helicase involved in DNA damage response and repair and adaptive immunity. It is an important research target in autoimmune disease and is widely studied to elucidate cellular mechanisms a



Immunocytochemistry analysis of Ku70 (green) in CEM using Ku70 antibody, and DAPI (blue).



Western blot analysis of Ku70 in HeLa, A549, MCF-7, U251 lysates using Ku70 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Ku70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.