

LIG1 Rabbit mAb


ABZ3327

Product Overview


Name	LIG1 Rabbit mAb
Accession(Primary)	P18858
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (P18858)


Gene Name	LIG1
About Gene Symbol	




Immunocytochemistry analysis of LIG1 (green) in Hela using LIG1 antibody, and DAPI(blue).



Western blot analysis of LIG1 in PC-12 lysates using LIG1 antibody.



Western blot analysis of LIG1 in HeLa, A549, HL-60 lysates using LIG1 antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using LIG1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.