

LILRB1 Rabbit mAb


ABZ3328

Product Overview

Name	LILRB1 Rabbit mAb
Accession(Primary)	Q8NHL6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 71 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:20-1:50; FC-1:200-1:500; IP-1:20-1:50

Target information (Q8NHL6)

Gene Name	LILRB1
About Gene Symbol	



Western blot analysis of LILRB1 in HL-60, U251 lysates using LILRB1 antibody.