

LINE 1 ORF1p Rabbit mAb


ABZ3329

Product Overview

Name	LINE 1 ORF1p Rabbit mAb
Accession(Primary)	Q9UN81
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (Q9UN81)

Gene Name	L1RE1
About Gene Symbol	



Western blot analysis of LINE1 ORF1p in HeLa, A549, HL-60, U2OS lysates using LINE 1 ORF1p antibody.