


LTK Rabbit mAb

ABZ3331


Product Overview

Name	LTK Rabbit mAb
Accession(Primary)	P29376
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 92 kDa; Observed MW: 92 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (P29376)

Gene Name LTK**About Gene Symbol**

Western blot analysis of LTK in A549, HL-60, U2OS, C6 lysates using LTK antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using LTK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.