

MAP3K4 Rabbit mAb


ABZ3332

Product Overview


Name	MAP3K4 Rabbit mAb
Accession(Primary)	Q9Y6R4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 182 kDa; Observed MW: 200 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:50-1:100

Target information (Q9Y6R4)


Gene Name	MAP3K4
About Gene Symbol	



Immunocytochemistry analysis of MAP3K4 (green) in Hela using MAP3K4 antibody, and DAPI (blue).



Western blot analysis of MAP3K4 in HL-60, U251 lysates using MAP3K4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using MAP3K4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.