

ME2 Rabbit mAb


ABZ3333

Product Overview


Name	ME2 Rabbit mAb
Accession(Primary)	P23368
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (P23368)


Gene Name	ME2
About Gene Symbol	



Immunocytochemistry analysis of ME2 (green) in Jurkat using ME2 antibody, and DAPI (blue).



Western blot analysis of ME2 in mouse intestines, rat spleen lysates using ME2 antibody.



Western blot analysis of ME2 in Ramos lysates
using ME2 antibody