

MELK Rabbit mAb


ABZ3334

Product Overview

Name	MELK Rabbit mAb
Accession(Primary)	Q14680
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 75 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q14680)

Gene Name	MELK
About Gene Symbol	MELK (Maternal Embryonic Leucine Zipper Kinase) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MELK in U2OS, C6 lysates using MELK antibody.