


MPP1 Rabbit mAb

ABZ3336

Product Overview

Name	MPP1 Rabbit mAb
Accession(Primary)	Q00013
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 52 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q00013)

Gene Name MPP1**About Gene Symbol**

Western blot analysis of MPP1 in A549, HL-60 lysates using MPP1 antibody.

