

MRAS Rabbit mAb


ABZ3337

Product Overview


Name	MRAS Rabbit mAb
Accession(Primary)	O14807
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (O14807)

Gene Name	MRAS
About Gene Symbol	MRAS (Muscle RAS Oncogene Homolog) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Mras in mouse brain lysates using MRAS antibody.



Western blot analysis of Mras in HeLa, rat Brain lysates using Mras antibody.