

NDUFS2 Rabbit mAb


ABZ3339

Product Overview


Name	NDUFS2 Rabbit mAb
Accession(Primary)	O75306
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:2000; IP-1:20-1:50

Target information (O75306)

Gene Name	NDUFS2
About Gene Symbol	



Western blot analysis of NDUF52 in HeLa, HepG2, U2OS, U251, HL-60 lysates using NDUF52 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using NDUF52 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.