

## NM23A Rabbit mAb


ABZ3340

### Product Overview

|                             |                                                       |
|-----------------------------|-------------------------------------------------------|
| <b>Name</b>                 | NM23A Rabbit mAb                                      |
| <b>Accession(Primary)</b>   | P15531                                                |
| <b>Host Species</b>         | Rabbit                                                |
| <b>Clonality</b>            | Monoclonal                                            |
| <b>Isotype</b>              | IgG                                                   |
| <b>Species Reactivity</b>   | Human                                                 |
| <b>Calculated MW</b>        | Calculated MW: 17 kDa; Observed MW: 17 kDa            |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>FC-1:50-1:100 |

### Target information (P15531)

|                          |                                                                                                                                                                                                                                                                                     |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | NME1                                                                                                                                                                                                                                                                                |
| <b>About Gene Symbol</b> | NME1 (NME/NM23 Nucleoside Diphosphate Kinase 1) is a protein kinase and signaling regulator involved in cell migration and invasion. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of NM23A in HeLa, A549, MCF-7, 293T lysates using NM23A antibody.