

NME2 Rabbit mAb

ABZ3341


Product Overview

Name	NME2 Rabbit mAb
Accession(Primary)	P22392
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; FC-1:20-1:50


Target information (P22392)

Gene Name	NME2
About Gene Symbol	


Mouse spleen



Western blot analysis of NME2 in mouse spleen lysates using NME2 antibody.



Western blot analysis of NME2 in U2OS, C6, Ramos lysates using NME2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using NME2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.