

NR2F2 Rabbit mAb


ABZ3342

Product Overview


Name	NR2F2 Rabbit mAb
Accession(Primary)	P24468
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P24468)

Gene Name	NR2F2
About Gene Symbol	NR2F2 (Nuclear Receptor Subfamily 2 Group F Member 2) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms.



Western blot analysis of NR2F2 in HeLa, A549, HL-60, U251, U87-MG lysates using NR2F2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using NR2F2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.